


Calorimeter EDM

by Denis Oliveira Damazio
(BNL)

Calorimeters

- * Different calorimeter technologies.
- * lead/Liquid Argon (EM - up to $\eta < 3.2$), copper/Liquid Argon (HEC - $1.5 < \eta < 3.1$) and copper-tungsten/Liquid Argon (FCAL - $3.0 < \eta < 4.9$).
- * Iron/Scintillating tiles (Hadronic - $\eta < 1.75$).
- * Highly segmented, specially in central EM region ($> 187k$).
- * Minimum Bias Trigger Scintillators (MBTS) using Tile hardware.
- * Gains to increase dynamic range : LAr (high, medium, low), Tile (high, low).
- * Used also for Trigger purposes (L1 and HLT).
- * Tile has two channels used for a single cell (two sides of the scintillator plates).
- * Upgrades trying to re-use a lot of the already developed code, specially EDMs.



How these calorimeter work ?

Check : <https://www.youtube.com/watch?v=WAHTMG5nMK4>

https://www.youtube.com/watch?v=mwpY7GamN_w

<https://www.youtube.com/watch?v=sj5QYCeQOOU>

Main classes



- * hardware-related classes : Raw Channels / Digits.
- * Physics related classes : CaloCell, LArCell, TileCell, geometry.
- * CaloCluster. Signal State change (look at different calibrations).

Digits

- * LAr and Tile Digits represent the sampled signal from the analogue calorimeter pulse.
- * usually only available in special MC productions or specific data samples, but also available in RAW for every cell above a given energy threshold.
- * Associates to hardware identifier and Gain and samples.

```
class LArDigit {
private:
    /** @brief Online Identifier */
    HWIdentifier m_hwIdentifier;
    /** @brief gain */
    CaloGain::CaloGain m_gain;
    /** @brief vector of ADC samples */
    std::vector<short> m_samples;
```

TileDigits :
hardware ID (with Gain information associated)
samples



Raw Channels

- * Representation of the hardware channel associated to a cell (two/cell for Tile) : Direct dump of Digital Signal Processing (DSP) units. No fancy calibration.
- * hardware ID connected to front-end channel (readout). Energy directly related to pulse amplitude, time related to shift (LHC jitter, time of flight), quality of the pulse fitting by optimal filter algorithm.
- * Usually available in RDOs (special production for trigger purposes).

LArRawChannel :

- * hardware ID (with Gain information associated)
- * energy
- * time (shift from reference)
- * quality/provenance (related to how well signal was fitted)
- * gain

TileRawChannel :

- * hardware ID (with Gain information associated)
- * amplitude
- * time (shift from reference)
- * quality (related to how well signal was fitted)
- * pedestal

CaloCell, LAr and Tile

- * Associates energy, time, quality to physics relevant position (η, ϕ , layer, detector type coordinates).
- * all coordinates come via geometry description (see next slide).
- * Many possible improvements on top of RawChannels (average energy for missing cells, pedestal train-bunch structure compensation, High-Voltage).
- * Full container available ("AllCalo") in ESDs level : no xAOD direct (ROOT-like) access.
- * Derived classes LArCell and TileCell used for very specific functions (e.g: access to independent Tile readouts). In Tile, "Gain" comes independently from each photomultipliers.

CaloCell :


- * Offline ID.
- * energy
- * time
- * quality/provenance
- * gain
- * CaloDetDescrElement
(geometry).

Geometry

- * Detector geometry has a raw version defined by xml dictionaries.
- * cells are later re-aligned from DB in the “Start of Run” event.
- * Precise position (η, ϕ , raw η , raw ϕ , x, y, z, r , raw x , raw y , raw z , raw r) and size ($\Delta\eta, \Delta\phi$, volume, $\Delta r, \Delta x, \Delta y, \Delta z$). Some methods (e.g: sampling layer) comes from yet another class (CaloDetDescriptor).
- * Because loading these conditions depends on run number->IOV access, cells access cannot happen at ROOT level.

Containers

- * LArRawChannels are organised by Front-End boards (128 channels) in LArRawChannelCollections. A Final LArRawChannelContainer contains all collections.
- * TileRawChannels are organised by Tile Drawer (46 channels - a phi-wedge) into TileRawChannelCollections. The TileRawChannelContainer contains all collections as well.
- * CaloCellContainer has all the cells (usually in ESD) and has fast cell access methods, specially if the container is complete (organisation by hash ID).




CaloCluster

- * CaloCluster available in xAOD. SignalState (calibrated, uncalibrated, etc) is available.
- * CaloCluster has global properties (η, ϕ , energy, mass, time) as well as quantities by layer (eSample, etaSample, phiSample). Also FourMon_t.
- * Multiple “size”s are defined to help differentiate between cluster types (like Sliding-window - window size, Topo Cluster - cluster growing thresholds).
- * possibility to loop over cells composing the cluster (no ROOT-access!).
- * Large number of Moments (variables associated to the cluster : cluster centroid, fraction of energy in Core, calibration, etc).


Calorimeter Reconstruction

- * The digits are used to recover the information about amplitude.
- * Detailed information about the detectors (obtained with special calibration data taking) is stored in DBs organised by hardware channel ID.
- * Many special classes defined for DB access.



Calorimeter Reconstruction

- * The acquisition chain starting from the front-end digitisation of the analog pulses.
- * DSPs produce RawChannels (and for cells above an energy threshold, Digits are attached).
- * Reconstruction produces CaloCells from RawChannels (or from Digits refitted at this stage). Many cell corrections are applied.
- * Clusters are built from cells and recorded into xAOD. Group of cells associated to each CaloCluster can be stored (but only accessed within athena!).




Clustering

- * two main approaches : Sliding window and topological clustering. Special calibration routines to each case.
- * Sliding window clusters are found when a pre-defined size window of cells (say 3x3) in the full detector (or only EM) has more than a threshold of energy and is a local maximum.
- * The window slides over the full detector.
- * Single layer (or groups of layers) can be used to seed (pre-cluster window finding). Later other layers can be added. Cluster position is energy weighted position of the cells in its window.

topological clustering

- * clusters do not have a fixed form.
- * initial (a) seeds are found when cell energy $> 4\sigma$ of noise (including pile-up noise).
- * cluster grows (b) as touching cells have energy $> 2\sigma$.
- * Final layer (c) added with cluster energy $> 0\sigma$.
- * Cluster positions derived from cells (usually energy weighted sums of cells properties- UpdateKine).
- * Special algorithms separate contributions (weights) of cells to two different clusters.
- * Adjustable thresholds. Resilience to pile-up.



Associations

- * Using online helpers for cells “hardware features” :

```
const HWIdentifier chid = (*cont_it)->channelID();  
HWIdentifier FebID = m_onlineHelper->feb_Id(chid); // LArOnlineID
```

- * Using offline helpers for cells “physics coordinates” or hash organisation. :

```
Identifier theCellID = aCell->ID();  
IdentifierHash theCellHashID = m_calo_id->calo_cell_hash(theCellID); //CaloCell_ID
```

- * Cabling (connect online world - rawchannel/digit - to offline - cells) :

```
const Identifier id = m_larCablingSvc->cnvToldentifier(chid); //LArCablingService  
HWIdentifier onlId=m_larcablingSvc->createSignalChannelID(id);
```

Change State

- * It is possible to change the state of a clusters and later come back to the original state :

```
{ // start scope
    CaloClusterChangeSignalStateList list;
    for(xAOD::CaloClusterContainer::const_iterator itr = m_pCaloClusterContainer->begin();
        itr!=m_pCaloClusterContainer->end(); ++itr)
        list.add( (*itr), xAOD::CaloCluster::UNCALIBRATED );
    std::cout << "CLUSTER ISSUE check the first now it should be uncalibrated : "
    << (*m_pCaloClusterContainer->begin())->et() << " " << (*(m_pCaloClusterContainer->begin()))->et()
    << std::endl;
} // destroy scope // cluster should come back to its original state of mind...
```

Calorimeter EDM

- * Calorimeter digits, rawchannels, cells, cluster presented in the order they appear in reconstruction.
- * More detailed information (like cells geometry) stored in DB. No ROOT access.
- * Clusters can be formed by different algorithms (sliding window or topological clustering) and accessed in xAOD level.

Event displays

